

Quickstart Guide

MilliKey™ LH-4 Button Box



Default MilliKey LH-4 Button to Keyboard Key Mappings

Thank you for choosing MilliKey.

The MilliKey LH-4 is a 1000 Hz USB Button Box that generates standard keyboard events when buttons are pressed.

Initial Set Up

■ Keyboard

The MilliKey is automatically detected as a 1000 Hz USB HID keyboard when connected to a Windows, macOS, or Linux computer. Start collecting standard keyboard data, it is that simple.

■ Serial

MilliKey button boxes also optionally utilize a 1000 Hz USB Serial interface, which is plug-and-play on Windows 10, macOS, and Linux. If using Windows 7 please install the required USB Serial driver from:

www.labhackers.com/downloads.html



LabHackers Inc.

www.labhackers.com

Quickstart Guide cont'd


MilliKey™ LH-4 Button Box

Modifying Settings

The MilliKey Manager application makes it easy to view and modify MilliKey settings such as button to keyboard key mappings.

Download MilliKey Manager for Windows from:

www.labhackers.com/downloads.html



Further Information

For any questions about the MilliKey LH-4, please contact:

Sales: sales@labhackers.com

Support: support@labhackers.com

For more information about other MilliKey models and options, please visit www.labhackers.com.



LabHackers Inc.

www.labhackers.com